

Place your logo here	Standard Operating Procedures
	Acknowledging Vital Sign Data

Purpose: This procedure provides a guideline for Assessment and Acknowledgement of Vital Sign Data in the Honeywell Life Care Solutions Monitoring System

Cross-Referencing: LifeStream User Manual, LifeStream Basics/The Current Status Window/Understanding the Current Status Area/ Understanding the Alert Condition Column/ Understanding Color Coded Alert Condition Indicators

Applications and Restrictions: Oversight of Vital Signs must be supervised by Clinical Management

Procedure Description:

1. Log into LifeStream
2. Double click on selected client for Tabular Trend Screen
3. Double click on line of Vital Sign data to open Respond to Reading Screen
 - A. Blue No Alert Parameters Set
 - Refer to Setting Alert Parameters
 - Respond to color coded data as indicated
 - B. Green WNL- all Vital Sign Data is within Alert Parameters set and the client answered 'No' to programmed question/s
 - Clinical assessment of trended data
 - Click on responses as indicated
 - Click on Acknowledge
 - C. Red Alert- a Vital Sign is outside of Alert Parameter set and/or answered Yes to programmed question/s
 - Clinical Assessment of trended data
 - Contact client/caregiver for additional subjective data
 - Determine appropriate interventions, request retest of Vital Signs, Contact Health Care Provider as per organization protocol or guideline
 - Document in 'Respond to Reading' Comment Section or click appropriate Responses
 - Click on 'Acknowledge' to complete or 'Hold' if indicated for additional follow up
 - D. Yellow Missing Data- indicates client data is missing for one or more vitals OR partial set of vitals received with retesting
 - Clinical assessment of trended data
 - Verify reason for missing data
 - Contact client if indicated to instruct them
 - Document in 'Respond to Reading' clicking appropriate Responses and free text in Comment Section as needed
 - Click on 'Acknowledge' to complete or 'Hold' if indicated for additional follow up
 - E. Orange Empty Packet-indicates monitor prompted client to take VS without response from client
 - Verify if client performed vital signs prior/post programmed time

- Verify on documentation that the client is unavailable due to physician appointment or for any other reason.
- Relative to findings contact client/contact person for follow up as per organization protocol or guideline
- Document in 'Respond to Reading' clicking appropriate Responses and free text in Comment Section as needed
- Click on 'Acknowledge' to complete or 'Hold' if indicated for additional follow up

** Test message sent – verify to installer that test message received on new installation

F. Orange Transmit Error-indicates no communication from monitor

- Verify possible reason for error
 - Monitor unplugged
 - Electricity off at programmed time
 - Monitor serial number may be different than the one assigned to client.
 - Incorrect PIN number in monitor
- Contact client/caregiver for information
- Troubleshoot reason for communication failure
- Contact Honeywell Customer Service as indicated
- Document in 'Respond to Reading' clicking on appropriate response and free text in Comment Section as needed
- Click on 'Acknowledge' to complete or 'Hold' if indicated for additional follow up

G. Invalidation of Partial or All Vital Sign Data

- For partial invalidation click on individual vital sign to uncheck
- Free text reason for invalidation in Comment Section
- Click on 'Save Validity' to keep remaining Vital Sign data valid
- Respond to remaining valid vital sign data as outlined above
- For invalidation of all Vital Sign Data free text reason in Comment Section and click on save.
- **Note:** Alert must be invalidated prior to being acknowledged

Steps for Maximizing Telehealth Alert Management:

1. **Establish your current total number of alerts.**
 - Establish a goal to decrease alerts from this baseline.
 - If you are at less than 50% per day, look at a realistic goal based on your client acuity and goal of your program.
2. **Ensure that staff installing the monitors follow best practices**
 - This includes sending a test message to ensure successful data transmission. ***As you know, not following this practice can delay client intervention(s) and creates empty packets.***
3. **Programmed questions should not exceed 3 to 4 on any given day per client.**
 - This assists in managing alerts and decreases client complacency when answering questions.
 - Consider programming different questions for a specific day(s) of the week. ***For example: "Have you been to the emergency room this week?" Ask on a Wednesday.***
4. **Explore changing any standard alert limit parameters:**
 - Some agencies will have their Professional Advisory Council or Medical Director Review alert limit parameters.
 - Also, ensure the policy/protocol indicates that if a client's normal parameters fall below the agency's standard, the telehealth clinician can alter the alerts based on the trends.
5. **Utilize trended data to adjust alert parameters at least weekly.** (This may require a physician's order.)
6. **Consider the use of smart limits.**
 - Agencies have utilized this feature as a starting point for setting alert parameters.