


LifeStream™ Web

What's New?

Cost-effectively and efficiently manage patient populations with a range of risk factors via a consolidated clinical dashboard.

Access LifeStream™ from the Web with Microsoft Internet Explorer®, version 11; or an updated, stable version, of the Google Chrome™ Web browser.




Video Visits

An intuitive interface makes it easy for LifeStream users to visit with patients, care advocates and family members. Leverage LifeStream to manage Video Visits from start to finish, invite participants and share content.




Custom Questions and Reminders

Add customized questions and reminders in your organization's database. Administrative edits to custom questions and reminders automatically apply to all monitor configurations and Presets.



Interactive Graphical Trends

Display an enhanced graphical trend report for any selected vital sign reading from the Dashboard, while showing all alert thresholds. Hover your mouse over any point in the graph to display reading details.



Reporting Capabilities

Query common reports from the Web, and print or export to common Adobe® Acrobat Reader® and Microsoft Excel® file formats for data sharing and analysis. New reports targeting hospitalization, hospitalization rate and monitor utilization data strengthen at-a-glance organizational data insights.



Advanced Scheduling Interface

A calendar-based interface helps users set and modify question schedules for compatible patient monitors. Browse a patient list by diagnosis, and deploy saved monitoring schedules to multiple patients from the same menu.

Technical Requirements

SOFTWARE REQUIREMENTS

Web Browser

Microsoft Internet Explorer® upgraded to Version 11

-or-

Google Chrome™ upgraded to the most recent, stable version.

Compatibility with other Web browsers has not been validated. Use of other web browsers may result in difficult-to-access controls or a cluttered display that may adversely impact ease of use.

Portable Document Format (PDF) Reader:

Use Adobe® Acrobat® Reader, or an updated PDF reader, to view downloaded files or user manuals.

Anti-Virus Software:

Install and maintain current anti-virus software on computers interacting with LifeStream™ to help protect patient data and your system.

Operating System:

LifeStream Web has been validated with a desktop Web browser, on a Microsoft Windows computer. Install and maintain all Operating System updates on any computers interacting with LifeStream™ to help protect patient data and your system.

SYSTEM REQUIREMENTS

Display: A monitor screen capable of supporting at least a 1200 x 768 resolution.

Peripherals: Mouse and keyboard.

Internet: An Internet connection is required.

SYSTEM RECOMMENDATIONS

Processor: 1GHz or faster x86- or x64-bit

Memory: 1GB RAM 32-bit; 2GB RAM 64-bit

Hard Drive: At least 1GB available

System Time: This web interface displays times in the time zone of the user's Web browser. To check the correct time for your time zone, use: <http://www.time.gov>.

LifeStream™ Updates: LifeStream™ uses a platform approach that limits the likelihood that incremental updates to either a supported operating system

or browser would impact the product. However, if a deleterious impact were reported, LifeStream™ supports an update mechanism that creates a path to deliver revised software. These impacts are limited to the specific client workstation and do not impact the retrospective collection of patient data.

Video Visits: A license is required to host video visits with a patient. All visits must include at least one patient. Genesis Touch® patient monitors must be upgraded to a compatible version, connected to Wi-Fi®, or licensed with a video-over-cellular data plan. Participants and hosts using Internet Explorer will be asked to download and install a Vidyo.io™ plugin, or a Vidyo.io™ application.

DEFINITIONS

Workstation: Any location where a user will access and interact with LifeStream™. Each workstation, both within the agency and those at any remote location, should be evaluated by your IT support staff to ensure the specifications have been met. The LifeStream™ Web interface may be accessed without an installation on a workstation, via a web interface.

Internet Access: The LifeStream™ Web client uses https:// connections. If an exception needs to be made in your Proxy Server or firewall to allow access, we recommend it be made to the domain vitalsmanager.com. Some network proxy configurations require authentication be entered at the web browser level, and may require that an exception be entered in the proxy for the system to work.

Some network proxy configurations require authentication be entered at the web browser level, and may require that an exception be entered in the proxy for the system to work.

HIPAA Medical Privacy Standards: LifeStream™ has been designed to comply with HIPAA regulations. It is the organization's responsibility for complying with privacy and security statutes for full HIPAA compliance, including any necessary and related policies and procedures. All users of the system hold responsibility and accountability for protecting the privacy of a client's health information.

Credentials: To install or access LifeStream™, authorized users need a username and password. To access the LifeStream™ Web interface, you may login at the web address given to you in your welcome documentation, or by your administrator.

Monitor Compatibility: LifeStream™ is compatible with Honeywell monitors and their compatible peripherals, except for Sentry II and Sentry 4-port pre-v2.04. Additionally, any interface product, SHP, or Vitals Manager currently working with Central Station v3.5 will no longer work with LifeStream™.

Export and Print: Install all printer drivers before attempting to print from LifeStream™.

Current Software: Maintain security and software updates provided by the manufacturer of your software. Certain software may no longer be supported. Microsoft® Windows XP® and Windows 2000® operating systems are no longer supported by Microsoft®. These operating systems can pose a cyber-risk via virus and Trojan horse infections that may exploit its weak-points.

HomMed®, LifeStream™ and Genesis Touch® are either trademarks or registered trademarks of Honeywell HomMed®, LLC. All other product names, trademarks, logos and brands cited herein, are property of their respective owners.

Android™, Google, Chrome™ browser and the Google logo is a trademark of Google Inc.; Adobe® and Adobe® Acrobat® Reader are registered trademarks of Adobe Systems Incorporated; Microsoft, Windows, Excel, Internet Explorer, Windows XP, Windows 7 and Windows 2000 are either registered trademarks or trademarks of Microsoft Corporation. Wi-Fi is a registered trademark of Wi-Fi® Alliance. VIDYO, the Vidyo cube icon, the Vidyo Partner Logos and names of Vidyo Inc.'s products and offerings are trademarks or registered trademarks of Vidyo, Inc., in the United States and in other jurisdictions.

Customer Service: 1-888-353-5404 | Ext. 1

Honeywell HomMed LLC is not responsible for the technical support of any workstations defined herein.

For more information

www.honeywelllifecare.com

Honeywell Life Care Solutions

3400 Intertech Drive, Suite 200
Brookfield, WI 53045 USA

E18-0004 | P2000EN.12B | 02/2018
© 2018 Honeywell HomMed LLC.

Honeywell